

Determination of polycyclic aromatic hydrocarbons and fatty acid methyl esters in diesel fuel by comprehensively coupled (supercritical fluid \times gas) chromatography

Daniel Malan

Department of Chemistry, University of Pretoria

Egmont Rohwer*

Department of Chemistry, University of Pretoria

Abstract

Historically, diesel fuel from petroleum sources consisted of aliphatic hydrocarbons, but it may now contain polycyclic aromatic hydrocarbon (PAH) and fatty acid methyl esters (FAME). PAHs might originate from the refining process, or might be blended in to increase the density or cetane number to comply with standards. FAMEs are blended in because it reduces the carbon footprint of the fuel, and because it helps meet the lubricity requirement.

Keywords: supercritical fluid chromatography, diesel, fatty acid methyl esters, polycyclic aromatic hydrocarbons, biodiesel

1. Introduction

Fuel for compression-ignition engines (diesel fuel) was traditionally obtained from petroleum sources. Today these traditional fuels are being blended with and replaced by diesel fuel from biological sources. To be sold on the open market, these fuels must comply with certain technical requirements, set by standards organizations such as ASTM and CEN. Because of

*Corresponding author

Email addresses: niel.malan@tuks.co.za (Daniel Malan),
egmont.rohwer@up.ac.za (Egmont Rohwer)

7 the complexity of diesel fuel obtained from petroleum sources (petrodiesel)
8 the chemical composition of the fuel is in general not prescribed. Instead,
9 a compliant fuel must pass certain distillation tests. There are two excep-
10 tions: the amount of polycyclic aromatic hydrocarbons (PAHs) and fatty
11 acid methyl esters (FAMES) are subject to upper limits.

12 PAHs are natural constituents of petrodiesel. They are noxious pollu-
13 tants, and therefore regulators try to reduce the exposure of humans to PAHs
14 by limiting their emission into the environment. Since a major source of PAH
15 emissions is internal combustion engines, one way to limit emissions is to limit
16 the amount in diesel fuels. A full discussion on PAH sources is beyond the
17 scope of this paper.

18 FAMES are not found in petroleum but are blended with petrodiesel.
19 FAMES are usually obtained by the transesterification of plant oils with
20 methanol. The resulting product can be used as a diesel fuel and is called
21 biodiesel. The primary benefit of the blending of FAMES is that it lowers
22 the carbon footprint of the fuel, but it can also improve the lubricity of high-
23 paraffinic diesel (1). Diesel for sale in the European Union can contain up to
24 7% FAMES. Any concentration of biodiesel makes a diesel fuel but, higher
25 concentrations are regulated as blends.

26 The amount of PAHs in diesel can be determined by a fluorescent indi-
27 cator absorption (FIA) test (ASTM D1319, (2)), but this method is labour-
28 intensive. Another method is a liquid chromatography method (EN 12916:2006)
29 with refractive index detection, and the method prescribed by ASTM D2425
30 requires a mass spectrometer. The method ASTM D5168 is a supercritical
31 fluid chromatography (SFC) method with flame ionization detection (FID).
32 Because the FID responds to number of carbon and has a large dynamic
33 range atoms in the sample FID methods are preferable for quantification.

34 The prescribed methods for determining FAMES in diesel is an infra-red
35 absorption method, ASTM D7806-12. This method is liable to interference
36 from FAMES degradation products that might accumulate during storage
37 (3).

38 The SFC-FID method ASTM 5186-19 claims that FAMES interfere with
39 the method. Weber et al. (4) claims that FAMES ‘foul’ the column. This
40 is contrary to experience: FAMES do elute from silica columns with neat
41 CO₂ as mobile phase (5). But the capacity factors (k' values) of the PAHs
42 are less than 0.05, while the capacity factors of the FAMES are greater than
43 0.8^{To do} (1). If a number of determinations of PAHs in FAME-containing
44 diesel fuel are done in sequence without allowing the FAMES to elute, the

45 FAMES will start eluting during the PAH determination and cause interfer-
46 ence. To address this problem The California Air Resources Board (ARB)
47 has published a standard operating procedure (SOP) (6) based on a publica-
48 tion by Diehl and DiSanzo (7) that combines the determination of PAHs with
49 FAMES. In this procedure SFC-FID is used to determine the PAHs, and then
50 the column is back-flushed to determine total FAMES. This solves the prob-
51 lem of FAMES interfering with the PAHs determination and simultaneously
52 quantifies the FAMES.

53 The information produced by this method is the total amount of FAMES
54 and the total amounts of mono-, di-, and tri+ aromatics. This information
55 suffices for monitoring compliance, but extending the method by comprehen-
56 sively coupling a fast gas chromatograph (GC) to the SFC provides a richly
57 patterned chromatogram that can help solve process control and supply chain
58 problems. The orthogonality of SFC and GC is high: on SFC with bare sil-
59 ica as a stationary phase and neat carbon dioxide as a mobile phase PAHs
60 are separated according to the number of rings in their structure, and the
61 FAMES are separated according to the number of double bonds, whereas all
62 GC separations are primarily according to volatility. Therefore determining
63 PAHs and FAMES make an ideal combination for comprehensively coupled
64 (supercritical fluid \times gas) chromatography. If — using the ARB SOP as a
65 basis — the FAMES are not back-flushed but all the compounds are allowed
66 to elute, and the SFC is comprehensively coupled to gas chromatography
67 (GC), the ¹D (SFC) dimension will separate the PAHs into their classes, and
68 the FAMES according to the number of double bonds. The ²D (GC) sepa-
69 ration will separate the compounds in each PAH and FAME class further,
70 according to their volatility. The resulting chromatogram contains patterns
71 that are rich in information.

72 **2. Material and methods**

73 **3. Results**

74 Figure 1 shows some results.

75 **4. Discussion**

76 **5. Conclusion**

77 **Acknowledgements**




Figure 1: . A contour plot of an SFCxGC chromatogram, showing the hydrocarbons in detail

78 **To do...**

79 □ 1 (p. 2): Use k' values from own experiments

80 **References**

- 81 [1] J. Rodríguez-Fernández, A. Ramos, J. Sánchez-Valdepeñas, J. R. Ser-
82 rano, Lubricity of paraffinic fuels additivated with conventional and
83 non-conventional methyl esters, *Advances in Mechanical Engineering* 11
84 (2019) 168781401987707.
- 85 [2] H. S. Knight, S. Groennings, Fluorescent indicator adsorption method for
86 hydrocarbon type analysis, *Analytical Chemistry* 28 (1956) 1949–1954.
- 87 [3] D. M. M. Pinho, V. O. Santos, V. M. L. dos Santos, M. C. S. Oliveira,
88 M. T. da Silva, P. G. T. Piza, A. C. Pinto, M. J. C. Rezende, P. A. Z.
89 Suarez, Evaluating the use of EN 14078 for determination of biodiesel in
90 diesel blends sold in the brazilian market, *Fuel* 136 (2014) 136–142.
- 91 [4] B. M. Weber, P. Walsh, J. J. Harynuk, Determination of hydrocarbon
92 group-type of diesel fuels by gas chromatography with vacuum ultraviolet
93 detection, *Analytical Chemistry* 88 (2016) 5809–5817.
- 94 [5] R. M. Smith, M. Hyytiäinen, A. V. Felipe, P. M. Morris, High resolution
95 packed column supercritical fluid chromatography of fatty acid methyl
96 esters, *Journal of Separation Science* 24 (2001) 208–212.
- 97 [6] Air Resources Board, Procedure for the analysis of biodiesel content in
98 diesel fuels by supercritical fluid chromatography and flame ionization
99 detector, 2017.
- 100 [7] J. W. Diehl, F. P. DiSanzo, Determination of total biodiesel fatty acid
101 methyl, ethyl esters, and hydrocarbon types in diesel fuels by supercritical
102 fluid chromatography-flame ionization detection, *Journal of Chromato-*
103 *graphic Science* 45 (2007) 690–693.